

The Pharmagard™ NTE 797 is a negative pressure total exhaust barrier isolator that meets the requirements of the **USP Std. 797** for sterile product preparation. **Meets: NIOSH & ASHP Hazardous Drug Preparation Stds.**

The NuAire Difference

- Class 100 [ISO level 5] clean air work area 99.99% HEPA filter
- True Laminar Airflow Minimizes air turbulence / cross contamination
- Excellent visibility into work area - expanded vision zone
- Two-part flexible work gloves
- Large glove ports allowing a larger range of motion - ergonomically enhanced
- Removable stainless steel work trays
- External cool white lighting reduces glare and interior heat build-up
- Large transfer chamber with sliding door
- Interior IV bar with two settings
- Mobile and adjustable with optional castored base stand
- Ability to sit or stand at a range of heights
- Will fit through standard door
- Fresh air enters cabinet through HEPA-filtered access port
- Exhaust air exits cabinet through HEPA-filtered exhaust port
- Requires building exhaust capability



Specifications:

Airflow Control System

Motor Speed Control:

Blower speed is regulated by solid-state motor voltage regulator

Supply HEPA filter:

Full size separatorless HEPA filter pack meeting the requirements of IEST-RP-CC-001 (type C)

Exhaust:

NU-NTE797-400: 215 CFM @ 1.0 w.g.
NU-NTE797-600: 305 CFM @ 1.0 w.g.

Installation:

Duplex Outlet:

Hospital Grade Duplex

Electrical Requirements:

115 VAC, 60 Hz, single phase, with a 10 amp maximum draw. Supplied with Hospital Grade line cord and plug. Listed by Underwriter's Laboratories

Front access:

Control panel, HEPA filters, and mini-helic pressure gauge are all front accessible.

Dimensions

Model NU-NTE797-400

Interior Dimensions:

(Height x Width x Depth)

Workzone:

27 1/2" x 36 1/4" x 24"

Transfer Chamber:

27 1/2" x 13 1/8" x 24"

Exterior Dimensions:

(Height x Width x Depth)

48" x 50" x 31"

Model NU-NTE797-600

Interior Dimensions:

(Height x Width x Depth)

Workzone:

27 1/2" x 60 1/4" x 24"

Transfer Chambers (2):

27 1/2" x 13 1/8" x 24"

Exterior Dimensions:

(Height x Width x Depth)

48" x 74" x 31"

Standard Features:

- HEPEX™ Zero Leak Airflow System
- Large Supply HEPA Filter
- Access and Exhaust HEPA Filters
- External Fluorescent Lighting
- Front Filter Removal
- Hinged 10 degree Slanted Window
- Removable Work Trays
- Minihelic Pressure Gauge
- Duplex Outlet
- Solid-State Motor Voltage Regulator
- Flexible Two Piece Glove Ports
- NU-797-400 - Two Glove Ports
- NU-797-600 - Three Glove Ports
- IV Bar With Two Settings
- Stainless Steel Interior

Optional Features and Accessories

- Ground Fault Interrupter
- Fixed Base Stand
- Motorized Adjustable Base Stand with 4" diameter castors
- Variety of Gloves
- Separate Sharps and Refuse disposal systems
- Gas Tight Damper
- Flex Duct Kits



Please note that performance specifications vary for models and types. Please consult NuAire for exact performance specifications

