

The Pharmagard™ PR797 is a positive pressure recirculating barrier isolator that meets the requirements of the USP Std. 797 for sterile product preparation.

### The NuAire Difference

- Class 100 [ISO level 5] clean air work area 99.99% HEPA filter
- True Laminar Airflow Minimizes air turbulence / cross contamination
- Excellent visibility into work area - expanded vision zone
- Two-part flexible work gloves
- Large glove ports allowing a larger range of motion - ergonomically enhanced
- Removable stainless steel work trays
- External cool white lighting reduces glare and interior heat build-up
- Large transfer chamber with sliding door
- Hinged access window offers complete access to work zone to facilitate easy cleaning
- Interior IV bar with two settings
- Mobile and adjustable with optional castored base stand
- Ability to sit or stand at a range of heights
- Will fit through standard door
- Fresh air enters cabinet through HEPA-filtered access port



### Specifications:

#### Airflow Control System

##### Motor Speed Control:

Blower speed is regulated by solid-state motor voltage regulator

##### Supply HEPA filter:

Full size separatorless HEPA filter pack meeting the requirements of IEST-RP-CC-001 (type C)

#### Installation:

##### Duplex Outlet:

Hospital Grade Duplex

##### Electrical Requirements:

115 VAC, 60 Hz, single phase, with a 10 amp maximum draw. Supplied with Hospital Grade line cord and plug. Listed by Underwriter's Laboratories

##### Front access:

Control panel, HEPA filters, and mini-helic pressure gauge are all front accessible.

#### Dimensions

##### Model NU- PR797-400

##### Interior Dimensions:

(Height x Width x Depth)

##### Workzone:

27 1/2" x 36 1/4" x 24"

##### Transfer Chamber:

27 1/2" x 13 1/8" x 24"

##### Exterior Dimensions:

(Height x Width x Depth)

48" x 50" x 31"

##### Model NU- PR797-400

##### Interior Dimensions:

(Height x Width x Depth)

##### Workzone:

27 1/2" x 60 1/4" x 24"

##### Transfer Chambers (2):

27 1/2" x 13 1/8" x 24"

##### Exterior Dimensions:

(Height x Width x Depth)

48" x 74" x 31"

#### Standard Features:

HEPEX™ Zero Leak Airflow System  
 Large Supply HEPA Filter  
 Access and Exhaust HEPA Filters  
 External Fluorescent Lighting  
 Front Filter Removal  
 Hinged 10 degree Slanted Window  
 Removable Work Trays  
 Minihelic Pressure Gauge  
 Duplex Outlet  
 Solid-State Motor Voltage Regulator  
 Flexible Two Piece Glove Ports  
 NU-797-400 - Two Glove Ports  
 NU-797-600 - Three Glove Ports  
 IV Bar With Two Settings  
 Stainless Steel Interior

#### Optional Features and Accessories

Ground Fault Interrupter  
 Fixed Base Stand  
 Motorized Adjustable Base Stand  
 with 4" diameter castors  
 Variety of Gloves  
 Separate Sharps and Refuse  
 disposal systems  
 Canopy Exhaust Transition  
 Flex Duck Kits



Please note that performance specifications vary for models and types. Please consult NuAire for exact performance specifications