

MICRONIC | SAMPLE STORAGE | TUBES

TUBESTRIP-8

DOES YOUR (CLINICAL) TESTING REQUIRE THAT YOUR SAMPLES ARE IN A GROUP OF 8? IF SO, THEN MICRONIC OFFERS THE TUBESTRIP-8, WHICH IS A LOW COST SOLUTION FOR STORING UP TO 8 SAMPLES IN A LINKED ROW. THE TUBESTRIP-8 HAS AN EFFECTIVE VOLUME OF 1.00ML (GROSS VOLUME 1.40ML). THE MICRONIC RANGE OF CAPS PROVIDES FOR THE SECURE CLOSURE OF THE TUBESTRIP-8.

Micronic offers the Tubestrip-8 a non-coded strip of 8 in a linked row of 1.40ml tubes. The tubes are manufactured of the highest-purity polypropylene, which resists many organic solvents.

The gross volume of the Tubestrip-8 is 1.40ml, net working volume is 1.00ml when used with Plug Style TPE or EVA caps. The Tubestrip-8 has a inner U-bottom.

When used in conjunction with Micronic tube caps or bands, the Tubestrip-8 provides for the simple and high integrity storage of clinical or biological samples, even at -80°C.


Manufactured from the highest grade pure polypropylene, the Tubestrip-8 has the same features and benefits as the complete range of Micronic sample storage tubes.

In combination with the TPE Capband-8 in 10 different colours or the economical EVA Capband-8, the Tubestrip-8 offers a complete solution for storing samples.

Offering full compatibility with Micronic tube racks, the new Tubestrip-8 integrates easily into automated laboratory environments and processes.

Available features

- Sterile, gamma radiated
- Snap, locks into the wells in the rack to prevent sample loss from overturned racks
- Pre-capped, with TPE Caps or Screw Caps



INFORMATION

- Gross volume 1.40ml
- Net volume 1.00ml when used with TPE or EVA caps
- RNase, DNase and Pyrogene free
- Class 7 clean room production
- Automation friendly
- Resists high chemical solvents




MICRONIC | SAMPLE STORAGE | TUBES

SPECIFICATIONS

MATERIAL	HIGHEST PURITY POLYPROPYLENE
GROSS VOLUME	1.40ML
TEMPERATURE RANGE	-80°C TO 121°C
RESISTANCE	RESISTS HIGH CHEMICAL SOLVENTS
CODING	NON-CODED



TECHNICAL DRAWING



ORDERING INFORMATION

■ 1.40ml Tubestrip-8		
MP32067	1.40ml Non-coded tubestrip-8 U-bottom	box/10
■ 1.40ml Tubestrip-8 Snap		
MP32067L	1.40ml Non-coded tubestrip-8 U-bottom Snap	box/10

Please visit our website www.Micronic.com or www.integra-biosciences.de



INTEGRA Biosciences GmbH
 Ruhberg 4 • 35463 Fernwald
 Tel. 06 404 / 809 / 0
 Fax 06 404 / 809 / 251
info@integra-biosciences.de
www.integra-biosciences.de

