

MICRONIC | SAMPLE STORAGE | RACKS

EXTENDABLE STORAGE HOLDER

DO YOU WANT A COST EFFECTIVE WAY TO STORE ANY STORAGE RACKS INTO A CHEST FREEZER? IF SO, THEN MICRONIC OFFERS YOU A FLEXIBLE AND LOW COST SOLUTION. THE MODULAR SYSTEM FITS IN ALL STANDARD STORAGE FREEZERS, IN TURN PROVIDING AN EFFECTIVE METHOD FOR STORING YOUR VALUABLE SAMPLES.

The Micronic Extendable Storage holder provides a storage solution for samples into a chest freezer. Simply 'click' the individual extendable storage holders together to a tower of max. 15 holders.

One holder has a capacity of 1x Comorack-96 or 1x Roborack-96, using the 1.10ml and 1.40ml tubes. When using the space saving 0.50ml or 0.75ml tubes in a Lobarack-96, the holder fits 2x Lobaracks-96, which doubles your storage capacity.


Using the Extendable Storage holder system saves valuable freezer space.

For the secure tracking of your samples, all the tubes are available with a laser etched, 2D Tracker Data-Matrix or TraXis code on the bottom of the tube.

The 2D Tracker Data-Matrix or TraXis coded tubes can be decoded with one of our 2D code readers:

The Single Rack Scanner (MVS) to decode 96 tubes in 15 seconds, and the single tube scanner to simply decode one tube at a time.

For High throughput application we offer a CCD camera based system, decoding 96 tubes within just 1 second.



INFORMATION

- Cost effective solution
- Lightweight storage system
- Compatible with standard storage racks
- Easy to use modular system
- Class 7 clean room production
- RNase & DNase and Pyrogene free




MICRONIC | SAMPLE STORAGE | RACKS

SPECIFICATION

MATERIAL	ABS-PC
RESISTANCE	RESISTS HIGH CHEMICAL SOLVENTS
AUTOCLAVABLE	YES
CAPACITY	1X COMORACK-96 OR ROBORACK-96 PER HOLDER 2X LOBORACK-96 PER HOLDER



TECHNICAL DRAWING



ORDERING INFORMATION

■ Extendable Storage Holders		
MP34000	Extendable Storage Holders	box/10

Please visit our website www.Micronic.com or www.integra-biosciences.de

