

# EXPERIENCE INNOVATIVE SOLUTIONS

## BIOLOGICAL SAFETY CABINET


BASE STANDS, EXHAUST TRANSITIONS, AND ERGONOMIC ACCESSORIES FOR  
IMPROVED SAFETY CONDITIONS, PRODUCTIVITY, AND QUALITY & RELIABILITY



## BIOLOGICAL SAFETY CABINET STANDS

ALL STANDS ARE FOR USE WITH LABGARD MODELS:  
 NU-425, NU-427, NU-430, NU-437, AND NU-440

NU-400-112  
 Telescoping w/Leg Levelers  
 30"-36" 4'



CATALOG NUMBER DESCRIPTION

### BASE STANDS—ADJUSTABLE

NU-400-111	Telescoping w/Leg Levelers 30"-36" 3'
NU-400-112	Telescoping w/Leg Levelers 30"-36" 4'
NU-400-113	Telescoping w/Leg Levelers 30"-36" 6'
NU-400-146	Telescoping w/Leg Levelers 30"-36" 5'

NU-400-128  
 30" w/4" Casters 4'



NU-400-127	30" w/4" Casters 3'
NU-400-128	30" w/4" Casters 4'
NU-400-129	30" w/4" Casters 6'
NU-400-148	30" w/4" Casters 5'
NU-400-131	36" w/4" Casters 3'
NU-400-132	36" w/4" Casters 4'
NU-400-133	36" w/4" Casters 6'
NU-400-149	36" w/4" Casters 5'

### ALL SUPPORT STANDS ARE CONSTRUCTED OF...

- welded 2" cold-rolled tubular steel with baked powder coat paint finish
- 30" Support Stand refers to the interior height of the work tray after complete assembly. Actual height of stand is 26½"
- 36" Support Stand refers to the interior height of the work tray after complete assembly. Actual height of stand is 32½"
- Leg Levelers adjustment: ± ¼"
- Telescoping base has two height settings—30" or 36" (work tray height)
- Manual crank or motorized adjustable actual standard height range 26"-34"

\* Factory Assembled (Shipped ready to use)  
 # Stands require assembly

Custom Base Support Stands with different heights are available

NU-400-141  
 Motorized Adjustable 4'




NU-400-140	Motorized Adjustable w/Leg Levelers 3'
NU-400-141	Motorized Adjustable w/Leg Levelers 4'
NU-400-142	Motorized Adjustable w/Leg Levelers 6'
NU-400-143	Motorized Adjustable w/Leg Levelers 5'

NU-400-156  
 Motorized Adjustable 4'  
 w/casters



NU-400-155	Motorized Adjustable w/ Casters 3'
NU-400-156	Motorized Adjustable w/ Casters 4'
NU-400-158	Motorized Adjustable w/ Casters 5'
NU-400-157	Motorized Adjustable w/ Casters 6'

## BIOLOGICAL SAFETY CABINET STANDS



*NU-400-104  
30" w/Leg Levelers 4'*

CATALOG NUMBER	DESCRIPTION
NU-400-102	30" w/Leg Levelers 2'
NU-400-103	30" w/Leg Levelers 3'
NU-400-104	30" w/Leg Levelers 4'
NU-400-105	30" w/Leg Levelers 6'
NU-400-144	30" w/Leg Levelers 5'
NU-400-107	36" w/Leg Levelers 2'
NU-400-108	36" w/Leg Levelers 3'
NU-400-109	36" w/Leg Levelers 4'
NU-400-110	36" w/Leg Levelers 6'
NU-400-145	36" w/Leg Levelers 5'



*NU-400-115  
Adjustable, LH Crank  
w/Leg Levelers 4'*

NU-400-114	Adjustable, LH Crank w/Leg Levelers 3'
NU-400-115	Adjustable, LH Crank w/Leg Levelers 4'
NU-400-116	Adjustable, LH Crank w/Leg Levelers 6'
NU-400-147	Adjustable, LH Crank w/Leg Levelers 5'



*NU-400-119  
30" w/Leg Levelers 6'*

### BASE STAND STORAGE CABINET

NU-400-117	30" w/Leg Levelers 3'
NU-400-118	30" w/Leg Levelers 4'
NU-400-119	30" w/Leg Levelers 6'
NU-400-151	30" w/Leg Levelers 5'
NU-400-120	36" w/Leg Levelers 3'
NU-400-121	36" w/Leg Levelers 4'
NU-400-122	36" w/Leg Levelers 6'
NU-400-152	36" w/Leg Levelers 5'




*NU-400-124  
LH Crank w/  
Leg Levelers 4'*

### BASE STAND STORAGE CABINET—ADJUSTABLE

NU-400-123	LH Crank w/Leg Levelers 3'
NU-400-124	LH Crank w/Leg Levelers 4'
NU-400-125	LH Crank w/Leg Levelers 6'
NU-400-153	LH Crank w/Leg Levelers 5'

## EXHAUST TRANSITION OPTIONS

ALL TRANSITIONS ARE FOR USE WITH LABGARD  
MODELS: NU-425, NU-437, AND NU-440



**GAS-TIGHT /  
LOW-PROFILE**



**BUTTERFLY VALVES**



**THIMBLE**



**SILICONE  
CONNECTION**



**8" TO 10" ROUND**



**FLEX DUCT**



**STAINLESS STEEL  
TURNTABLE**



**EXHAUST AIR FLOW  
MONITOR**



**ADJUSTABLE  
FOOTREST**



**REMOTE-CONTROLLED  
SERVICE VALVES**



**ELBOW RESTS**



**FOAM ARM PAD**

## ERGONOMIC FEATURES AND ACCESSORIES

AUTOMATIC ADJUSTABLE BASE SUPPORT STAND • ADJUSTABLE FOOTREST • ELBOW REST  
FOAM ARMREST PADS • STAINLESS STEEL TURNTABLE • REMOTE-CONTROLLED SERVICE VALVES



[www.nuaire.com](http://www.nuaire.com)

2100 Fernbrook Lane  
Plymouth, MN 55447  
U.S.A.  
Fax: 763.553.0459  
1.800.328.3352